

ddpcRquant Hands-on

Wim Trypsteen

Matthijs Vynck


Doctoral schools specialist course

1 september 2016, 't Pand Ghent

Day 1 (01-09-2016)	Day 2 (02-09-2016)	
08:00-9:30	Introduction to digital PCR <i>Ward De Spiegelaere</i>	Assay Design + Hands on training <i>Jo Vandesompele</i>
09:30-10:00	coffee break	coffee break
10:00-12:00	Absolute quantification by dPCR: theory + applications (HIV, excision circles, splice variants) <i>Linos Vandekerckhove</i> <i>Ide Smets</i>	Multitarget applications: CNV, human and mouse examples, genotyping, mutation detection <i>Jo Vandesompele</i> <i>Nicole Hersmus</i> <i>Patrick Pauwels</i>
12:00	Lunch	Lunch
13:00-14:30	Statistical considerations of dPCR and droplet calling methods <i>Lieven Clement</i> <i>Olivier Thas</i>	RNA quantification by dPCR and normalization tools <i>Ward De Spiegelaere</i> <i>Nicole Hersmus</i>
14:30-15:00	coffee break	coffee break
15:00-17:00	exercices on Shiny app <i>Wim Trypsteen en Matthijs Vynck</i>	GLMM models for dPCR + Hands on <i>Lieven Clement</i> <i>Olivier Thas</i>
17:00-17:15	short break	Short break
17:15-18:00	Power analysis for dPCR <i>Olivier Thas</i>	Tips and Tricks Minimal reporting guidelines <i>Ward De Spiegelaere</i>

GOALS

1. Introduce everyone to the ddpcRquant tool
2. Independent analysis of ddPCR data



ddpcRquant overview

Annotation Requirements

Demo + DIY Analysis

Exporting Data from the ddpcr Reader

Running ddpcRquant

Exercises part 1 and 2

Head file options

Step1: Input & Automatization

Quantasoft (Bio-Rad)

Head file

Amplitude files

Well

Sample name

Type Assay

Assay

Identify wells for NTC and Samples per assay

Assay

NTC Wells

Sample Wells

Individual annotated wells

Step2: NTC pre-processing

Baseline correction

Merge multiple NTCs

Step3: NTC processing

Apply extreme value theory

Calculate threshold

Step4: Samples processing

Baseline correction

Apply threshold

Calculate concentration

Step5: Output

Plots

NTC threshold plot

Sample dataplot + threshold

Summary file

Info: well, assay, sample

Pos/Neg/Tot Droplets

Concentration + LL/UL

HEAD file

- Summary of the experiment
- Contains annotation info (NTC,...)



Well	ExptType	Experiment	Sample	TypeAssay	Assay	Status
A01	Absolute Quantification	Absolute Quantitation template	dil 160	Ch1Unknown	RU5	OK
A02	Absolute Quantification	Absolute Quantitation template	dil 10	Ch1Unknown	RU5	OK
A03	Absolute Quantification	Absolute Quantitation template	dil 160	Ch1Unknown	LTR GAG PETRA	OK

AMPLITUDE files

- Individual well files with the fluorescent intensity information (droplets)



Assay1 Amplitude	Assay2 Amplitude	Cluster
1057.41455		1
1205.11		1
1227.16284		1
1266.01575		1
1290.73767		1

Step1: Input & Automatization

Quantasoft (Bio-Rad)

Head file

Amplitude files

Well

Sample name

Type Assay

Assay

Identify wells for NTC and Samples per assay

Assay

NTC Wells

Sample Wells

Individual annotated wells

Step2: NTC pre-processing

Baseline correction

Merge multiple NTCs

Step3: NTC processing

Apply extreme value theory

Calculate threshold

Step4: Samples processing

Baseline correction

Apply threshold

Calculate concentration

Step5: Output

Plots

NTC threshold plot


Sample dataplot + threshold

Summary file

Info: well, assay, sample

Pos/Neg/Tot Droplets

Concentration + LL/UL



Step1: Input & Automatization

Quantasoft (Bio-Rad)

Head file

Amplitude files

Well

Sample name

Type Assay

Assay

Identify wells for NTC and Samples per assay

Assay

NTC Wells

Sample Wells

Individual annotated wells

Step2: NTC pre-processing

Baseline correction

Merge multiple NTCs

Step3: NTC processing

Apply extreme value theory

Calculate threshold

Step4: Samples processing

Baseline correction

Apply threshold

Calculate concentration

Step5: Output

Plots

Summary file

NTC threshold plot

Info: well, assay, sample

Sample dataplot + threshold



Pos/Neg/Tot Droplets

Concentration + LL/UL

Fisher-Tippett theorem: The distribution of block maxima is given by the Generalized Extreme Value distribution (GEV)

= Block maxima (extremes) follow this family of distributions

merged_ntc_threshold_RU5.png



Step1: Input & Automatization

Quantasoft (Bio-Rad)

Head file

Amplitude files

Well

Sample name

Type Assay

Assay

Identify wells for NTC and Samples per assay

Assay

NTC Wells

Sample Wells

Individual annotated wells

Step2: NTC pre-processing

Baseline correction

Merge multiple NTCs

Step3: NTC processing

Apply extreme value theory

Calculate threshold

Step4: Samples processing

Baseline correction

Apply threshold

Calculate concentration

Step5: Output

Plots

Summary file





NTC threshold plot

Info: well, assay, sample

Sample dataplot + threshold

Pos/Neg/Tot Droplets

Concentration + LL/UL



$$C = -\ln\left(\frac{N_{neg}}{N}\right) * \frac{1000}{V_d} * D$$

D: Template/Mix

Well	assay	name	type	positive droplets	negative droplets	total droplets	concentration	lowerCI	upperCI
1	merged	2 LTR merged_NTC	ntc	1	19376	19377	0.5671	0.0666	4.8298
2	E01	2 LTR gDNA 1	sample	116	10957	11073	115.7274	91.1503	146.887
3	F01	2 LTR gDNA2	sample	163	10391	10554	171.0427	139.8199	209.1717
4	G01	2 LTR plas 1	sample	1	12425	12426	0.8844	0.1038	7.531
5	H01	2 LTR plas 2	sample	9	7758	7767	12.7409	5.5316	29.3334

ddpcRquant is made for 1-D analysis (Bio-Rad qx100/qx200)

- 1 channel at a time
- If 2 channels per well => one after the other will be analyzed

Use the NTCs to calculate a single threshold

- Whether the extreme value measured in the NTC is a real or false positive value is beyond the scope of ddpcrquant

ddpcRquant uses a predefined annotation

- Throws most of the errors when using the tool
- R functions available but webtool GUI

ddpcRquant overview

Annotation Requirements

Demo + DIY Analysis

Exporting Data from the ddpcr Reader

Running ddpcRquant

Exercises part 1 and 2

Head file options




Annotation before run starts: Quantasoft



Annotation after run

- Change headfile in text format reader
notepad, noteblock, NOT excel
- Use qx100/qx200 template



Before the ddPCR run: Annotation & Input Info

SampleName: NTC/ntc can only be in the name of the negative control

AssayName: different assay names required

Type: Select Ch1Unknown or Ch2Unknown (nothing else at the moment)

WellNumber: automatically



ddpcRquant overview

Annotation Requirements

Demo + DIY Analysis














Exporting Data from the ddpcr Reader

Running ddpcRquant

Exercises part 1 and 2

Head file options

After ddPCR run: Load ddPCR data into Quantasoft

		2013-03-29 plasmid dilutions 2LTR and Total DNA.bkp	29/03/2013 14:52	BKP-bestand	8594 kB
head ←		2013-03-29 plasmid dilutions 2LTR and Total DNA.csv	29/03/2013 14:52	CSV-bestand van ...	4 kB
		2013-03-29 plasmid dilutions 2LTR and Total DNA.log	29/03/2013 14:52	Tekstdocument	474 kB
QS ←		2013-03-29 plasmid dilutions 2LTR and Total DNA.qlp	29/03/2013 14:53	QuantaSoft Plate	8594 kB
		2013-03-29 plasmid dilutions 2LTR and Total DNA.qlr	29/03/2013 14:50	QLR-bestand	8547 kB
		2013-03-29 plasmid dilutions 2LTR and Total DNA_A01_RAW.qlb	29/03/2013 14:31	QLB-bestand	8786 kB
		2013-03-29 plasmid dilutions 2LTR and Total DNA_A02_RAW.qlb	29/03/2013 14:41	QLB-bestand	8770 kB
		2013-03-29 plasmid dilutions 2LTR and Total DNA_B01_RAW.qlb	29/03/2013 14:32	QLB-bestand	8770 kB
Raw data		2013-03-29 plasmid dilutions 2LTR and Total DNA_B02_RAW.qlb	29/03/2013 14:42	QLB-bestand	8770 kB
		2013-03-29 plasmid dilutions 2LTR and Total DNA_C01_RAW.qlb	29/03/2013 14:34	QLB-bestand	8770 kB
		2013-03-29 plasmid dilutions 2LTR and Total DNA_C02_RAW.qlb	29/03/2013 14:44	QLB-bestand	8770 kB
		2013-03-29 plasmid dilutions 2LTR and Total DNA_D01_RAW.qlb	29/03/2013 14:35	QLB-bestand	8770 kB
		2013-03-29 plasmid dilutions 2LTR and Total DNA_D02_RAW.qlb	29/03/2013 14:45	QLB-bestand	8770 kB

QuantaSoft Version 1.3.2.0

2013-09-18patients2LTRPicado Options

Hides the additional user options

BIO-RAD

Setup

Run

Analyze

About

Contact Support

Setup

Plate

Load Save As

Experiments

ABS
RED
CNV1
CNV2
CNV3
CNV4
Absolute Quantification

New Edit Remove

Template

New Load Save As

	01	02	03	04	05	06	07
A	pt1 A Absolute Quan... U 0.703	pt1 A Absolute Quan... U 0.307	pt1 A Absolute Quan... U 0.102	ntc gDNA A Absolute Quan... U 0			
B	pt2 A Absolute Quan... U 0.704	pt2 A Absolute Quan... U 0	pt2 A Absolute Quan... U 0.111	ntc gDNA A Absolute Quan... U 0			
C	pt3 A Absolute Quan... U 0.452	pt3 A Absolute Quan... U 0.169	pt3 A Absolute Quan... U 0.114	ntc gDNA A Absolute Quan... U 0			
D	pt4 A Absolute Quan... U 1.25	pt4 A Absolute Quan... U 0.0923	pt4 A Absolute Quan... U 0	ntc gDNA A Absolute Quan... U 0.0814			

QuantaSoft Version 1.3.2.0

2013-09-18patients2LTRPicado

BIO-RAD

Setup

Plate: Load, Save As

Template: New, Load, Save As

Experiments: ABS, RED, CNV1, CNV2, CNV3, CNV4, Absolute Quantification

Hide Options

Reprocess Raw Data

Export Amplitude and Cluster Data

Select Wells by Row

Charts

	01	02	03	04	05	06	07
A	pt1 Absolute Quan... 0.703	pt1 Absolute Quan... 0.307	pt1 Absolute Quan... 0.102	ntc gDNA Absolute Quan... 0			
B	pt2 Absolute Quan... 0.704	pt2 Absolute Quan... 0	pt2 Absolute Quan... 0.111	ntc gDNA Absolute Quan... 0			
C	pt3 Absolute Quan... 0.452	pt3 Absolute Quan... 0.169	pt3 Absolute Quan... 0.114	ntc gDNA Absolute Quan... 0			
D	pt4 Absolute Quan... 1.25	pt4 Absolute Quan... 0.0923	pt4 Absolute Quan... 0	ntc gDNA Absolute Quan... 0.0814			




















Run

Analyze

About

Contact Support

After exporting the folder should look like:

		2013-03-29 plasmid dilutions 2LTR and Total DNA.bkp	29/03/2013 14:52	BKP-bestand	8594 kB
head	←		2013-03-29 plasmid dilutions 2LTR and Total DNA.csv	29/03/2013 14:52	CSV-bestand van ... 4 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA.log	29/03/2013 14:52	Tekstdocument 474 kB
QS	←		2013-03-29 plasmid dilutions 2LTR and Total DNA.qlp	29/03/2013 14:53	QuantaSoft Plate 8594 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA.qlr	29/03/2013 14:50	QLR-bestand 8547 kB
Amplitude csv	←		2013-03-29 plasmid dilutions 2LTR and Total DNA_A01_Amplitude.csv	14/01/2014 15:12	CSV-bestand van ... 219 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA_A01_RAW.qlb	29/03/2013 14:31	QLB-bestand 8786 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA_A02_Amplitude.csv	14/01/2014 15:12	CSV-bestand van ... 165 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA_A02_RAW.qlb	29/03/2013 14:41	QLB-bestand 8770 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA_B01_Amplitude.csv	14/01/2014 15:12	CSV-bestand van ... 172 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA_B01_RAW.qlb	29/03/2013 14:32	QLB-bestand 8770 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA_B02_Amplitude.csv	14/01/2014 15:12	CSV-bestand van ... 163 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA_B02_RAW.qlb	29/03/2013 14:42	QLB-bestand 8770 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA_C01_Amplitude.csv	14/01/2014 15:12	CSV-bestand van ... 181 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA_C01_RAW.qlb	29/03/2013 14:34	QLB-bestand 8770 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA_C02_Amplitude.csv	14/01/2014 15:12	CSV-bestand van ... 141 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA_C02_RAW.qlb	29/03/2013 14:44	QLB-bestand 8770 kB
	←		2013-03-29 plasmid dilutions 2LTR and Total DNA_D01_Amplitude.csv	14/01/2014 15:12	CSV-bestand van ... 186 kB
			2013-03-29 plasmid dilutions 2LTR and Total DNA_D01_RAW.qlb	29/03/2013 14:35	QLB-bestand 8770 kB

This is how the csv file should look like

	A	B	C	D	E
1	Assay1 Amplitude, Assay2 Amplitude, Cluster				
2	1133.30322,,1				
3	1153.52869,,1				
4	1170.84326,,1				
5	1170.69263,,1				
6	1171.05542,,1				
7	1177.48816,,1				
8	1184.29419,,1				
9	1185.09314,,1				
10	1186.19812,,1				
11	1186.13586,,1				
12	1188.953,,1				
13	1189.22473,,1				
14	1190.09875,,1				
15	1192.60938,,1				
16	1192.57556,,1				
17	1193.63965,,1				
18	1194.11145,,1				
19	1194.71948,,1				
20	1195.99414,,1				
21	1195.99084,,1				
22	1196.14392,,1				
23	1198.23364,1				

2013-11-07 test_A01_Amplitude

On the first row:
3 names separated by 2 commas

Next rows:
3 numerical values separated by 2 commas

Depending on the measured channels 1 or 2 decimal numbers

Cluster number

Decimal sign needs to be “ . ”

Control Panel > Regional and Language options > Format tab

Land en taal

Notaties | Locatie | Toetsenborden en talen | Beheer

Notatie:
Nederlands (België)

Datum- en tijdsnotaties

Korte datumnotatie: d/MM/yyyy

Lange datumnotatie: dddd d MMMM yyyy

Korte tijdnotatie: U:mm

Lange tijdnotatie: U:mm:ss

Eerste dag van de week: maandag

[Wat houdt de notatie in?](#)

Voorbeelden

Korte datumnotatie:	1/06/2015
Lange datumnotatie:	maandag 1 juni 2015
Korte tijdnotatie:	10:28
Lange tijdnotatie:	10:28:49

[Meer instellingen...](#)

[Meer informatie online over het wijzigen van talen en regionale notaties](#)

OK Annuleren Toepassen

Indeling aanpassen

Getallen | Valuta | Tijd | Datum

Voorbeeld

Positief: 123456789.00 Negatief: -123456789.00

Decimaalteken: .

Aantal decimalen: 2

Cijfergroeperingssymbool: ,

Cijfergroepering: 123456789

Symbool voor negatief getal: -

Notatie van negatief getal: -1.1

Voorloopnullen weergeven: 0.7

Lijstscheidingstekens: ;

Maatstelsel: Metrisch

Standaardcijfers: 0123456789

Taalspecifieke cijfers gebruiken: Nooit

Klik op Beginwaarden als u de standaardinstellingen voor het systeem wat betreft getallen, valuta's, tijden en datums wilt instellen.

Beginwaarden

ddpcRquant overview

Annotation Requirements

Demo + DIY Analysis

Exporting Data from the ddpcr Reader

Running ddpcRquant

Exercises part 1 and 2

Head file options

Where to find the tool and additional info:

1. Browse to <http://www.ddpcrquant.ugent.be>
 - Click on the ddpcRquant webtool link (demo)
2. See also the Quick ddpcRquant analysis guide

What you will need today:

1. There is a preloaded demo dataset
2. Download the exercises folder from the course dropbox

ddpcRquant overview

Annotation Requirements

Demo + DIY Analysis

Exporting Data from the ddpcr Reader

Running ddpcRquant

Exercises part 1 and 2

Head file options

1. Perform the demo analysis (set to qx100!)
2. Use the demo dataset and set the threshold to a manual defined threshold for assay RNA
 1. Download the plots, summary table and summary parameters



1. Quantasoft analysis in following folder exercises_part2

1. Open in Quantasoft (qlp)
2. Go to Analyze (left side panel)
 1. Select all samples (*) + sort on target + select all RU5 + update
 2. Go to 1-D Amplitude
 3. Auto Analyze with individual and combined wells => What is wrong?
 1. Individual wells?
 2. Combined wells?
 4. Select wells C01 and E03 and look in ch1 => What is wrong?
 5. Write down the concentration of C01
3. Export the amplitudefiles
 1. Go to Setup > options > export amplitude and cluster data
 2. Select folder



2. Upload the amplitudefiles + headfile to ddpcRquant

3. Run analysis with standard settings

Write down the concentration of sample C01 (0.5 vs 0.28) => reason?

4. Run analysis with a 99.99 threshold

What happens with the threshold? Why?



ddpcRquant overview

Annotation Requirements

Demo + DIY Analysis

Exporting Data from the ddpcr Reader

Running ddpcRquant

Exercises part 1 and 2

Head file options



Annotation before run starts: Quantasoft



Annotation after run

- Change headfile in text format reader:
notepad, noteblock, NOT excel
- Use qx100/qx200 template

Open following headfile:



1. Find some missing or wrong annotation and change it according to the annotation rules

6 mistakes

2. Upload head & amplitude and look at overview

Demo: Make template



1. Make the demo template yourself

Wells: A03,B03,G03,H03

Sample: test1,test1,NTC,NTC




TypeAssay: Ch1Unknown

Assay: RPP30

2. Upload template & amplitudefiles to ddpcRquant

make sure the template file is uploaded first (sort folder before upload)

SIMPLY EXPLAINED



CONCURRENCY

Good luck with future analysis!

Troubleshooting + ddPCR advice:

Wim.trypsteen@ugent.be