

Quick Guide for ddpcRquant analysis

1. Exporting amplitudefiles from Quantasoft
 - Open the ddpcr main file in quantasoft: .qlp
 - Select the desired samples
 - Select options > export cluster and amplitude data
 - Press export

 - Make sure PC settings for a Decimal sign are set to a “ . ”
2. Upload the headfile and amplitudefiles to ddpcrquant
 - Plate annotation available: use the headfile (.csv)
 - No plate annotation available: add annotation after run
 - Change the headfile with notepad or other text based program (no excel)
 - Use the headfile template qx200/qx100 (excel)
 - Browse to <http://www.ddpcrquant.ugent.be> and click on the webtool link
 - Make sure your headfile is uploaded as the first file
 - Sort the folder based on name to get the headfile above the amplitudefiles
3. Go through the ddpcRquant workflow and download your data
 - Select your options in the NTC threshold analysis tab
 - Select Assay
 - Template volume and total mix volume: concentration calculations with dilution factor
 - Set both to 20 μ L (= dilution factor of 1) for concentration without dilution factor
 - Select threshold CI or give manual threshold
 - For each new dataset reset the analysis
4. Please cite ddpcRquant in your research articles!

Trypsteen et al. ddpcRquant: threshold determination for single channel droplet digital PCR experiments. Anal Bioanal Chem (2015) 407:5827–5834. PMID: 26022094